EN

Description

Adjustable glasses comprising 2 movable parts per lens. The power of each lens changes when the 2 parts move relative to one another.

Indications

This device is intended for long sightedness, short sightedness and presbyopia only. See specifications for power range.

Specifications

Lens power: -6.0 to +3.0 dioptres
Frame measurement: 49□22, 141 mm
Centration distance: 63 mm

Contents

Glasses with case: 1
Instructions for Use: 1

Directions for Use

Power Adjustment

1. Turn both dials forward and place glasses on head.
2. Cover left eye and turn right dial backwards until focus is achieved.
3. Make fine adjustments by turning dial forwards or backwards.
4. Cover right eye and turn left dial backwards until focus is achieved.
5. Make fine adjustments by turning dial forwards or backwards.
6. Your glasses are ready for use.

Disposal of devices

Dispose of device in accordance with local regulations.

Storage

Store this device in its protective case when not in use.

Cleaning Instructions

Use a lens cloth dampened with lens solution or soapy water to clean lenses. Dry lenses by rubbing lenses gently with a clean soft cloth (do not use tissue). Do not submerge lenses in liquid.

Warning

- This device is not intended to replace the need for regular examinations of eye health.
- This device should not be used for high-risk tasks.
- Not for use as eye protection.

Cautions

- Should this device cause any discomfort, do not wear.
- This device is not intended for astigmatism.
- This is a temporary device designed for short term use.
- Handle and clean with care.

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